Add-on IV Bag Spike with Integrated Clave®

For use with IV Admin set

EN	- DIRECTIONS
ZH	- 方向
DA	- BRUGSANVISNING
FI	- OJHEET
FR	- MODE D'EMPLOI
DE	- GEBRAUCHSANWEISUNG
IT	- ISTRUZIONI
JP	- 方向
NO	- ANVISNING
PT	- INSTRUÇÕES
ES	- INSTRUCCIONES
SE	- BRUGSANVISNING

US and Foreign Patents Pending

DF-3873 Rev 09 (03/18)

 ICU Medical, Inc.
951 Calle Amanecer
San Clemente, CA 92673 USA
Tel: +1 949 366 2183 Fax: +1 949 366 8368

EC REP
MDSS GmbH
Schiffgraben 41
30175 Hannover, Germany

STERILE|R

EN - English

Add-on IV Bag Spike with Integrated Clave®

For use with IV Admin set.

Not Made with DEHP.

Use Aseptic Technique. Single use only. Do Not resterilize.
Sterile, Non-pyrogenic fluid path in unopened undamaged package.

INDICATIONS FOR USE:

The ChemoClave Closed System Drug Transfer Device prevents the transfer of environmental contaminants, including bacterial and airborne contaminants into the system, and the escape of drug or vapor concentrations outside the system. The ChemoClave is needlefree and cannot be deactivated, which will passively aid in preventing needlestick injuries and the exposure to cytotoxic medications for healthcare personnel.

DIRECTIONS:

1. Prepare IV solution container in accordance with manufacturer's instructions.
2. Remove end cap and insert spike of pump admin set to the distal end of the ChemoClave IV Bag Spike. Close Clamps. If utilizing dry spike adapter: spike pump set straight into dry spike until membrane is punctured and rotate pump set with one full turn in one direction. Do not twist back and forth.
3. Remove end cap from bag spike and insert into solution container.
4. Suspend solution container and open clamp(s) to prime administration set. Repeat for multiple spikes if necessary.
5. Attach IV set to desired vascular access device for infusion.
6. (if applicable) Remove protective cap from Clave.
7. To use Clave spike port, swab and attach syringe to Clave by pushing and twisting $\frac{1}{4}$ turn until secure. Once connected, the Clave is open to a dedicated lumen for infusion or withdrawal directly with the solution container.
8. Replace protective cap on Clave with a cap that does not contain a male luer or cone inside the cap. (if applicable).

NOTES:

- A CSTD must not be used as a substitute for a containment primary engineering control (C-PEC). A CSTD may offer additional levels of protection and may facilitate enhanced occupational protection, especially when handling hazardous drugs outside of primary and secondary engineering controls (e.g., during administration).
- The ability to prevent microbial ingress should not be interpreted as modifying, extending, or superseding a drug manufacturer's labeling recommendations for the storage and expiration dating. Refer to the drug manufacturer's recommendations for shelf life and sterility information.

CAUTIONS:

- Tighten all connections before use.
- Do not use needles or luer caps on Clave.
- The Clave Connector is compatible with luers with an internal diameter (ID) between 1,55 mm and 2,8 mm.
- Use only ISO 8536-4:2010 compliant spikes
- Reuse negatively impacts performance / sterility potentially resulting in product failure / contamination.
- This device should be changed in accordance with current, recognized guidelines of IV therapy.
- Federal (USA) law restricts this device to sale by or on the order of a physician.

These directions for use were drawn up in English and have been translated into other languages. In the event of inconsistencies of interpretations, the English text shall be determinative.

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Directions for Use

Single patient use, do not resterilize. Sterile, non-pyrogenic fluid path in unopened undamaged package.

INDICATIONS FOR USE:

For use with vascular access devices, IV administration and extension sets for the direct injection, continuous or intermittent infusion, aspiration and needlefree delivery of drugs and solutions during IV therapy.

DIRECTIONS:

1. Using aseptic technique, open package and tighten connections where applicable.
2. For standalone Spiros, attach to male luer of administration set or syringe. Twist devices together until a click is heard and the Spiros begins to rotate freely. Spiros is now permanently attached to luer.
3. To prime and use various components, refer to chart below.
4. Once primed, attach to patient's catheter.
 - rotating luers, ensure the luer is fully engaged prior to locking down the collar;
 - locking blunt cannula connectors, connect to compatible pre-slit/pre-pierced injection port;
 - graduated adapters, push and twist until connector is secure.

NOTES:

- Pressure infusion tubing sets or removable extension sets made with pressure infusion tubing are noted in the name/description of the product. Pressure infusion tubing is rated to 400 psig. Pressure infusion tubing or extension sets should be connected to pressure infusion compatible IV sets only. For sets equipped with Clave/ MicroClave Y site: prior to pressurizing, activate the slide clamp and give the pressure infusion through the Clave/MicroClave Y site.
- Disinfect all needlefree ports using an aggressive circular motion for three (3) seconds.
- Do not leave open ports/hubs exposed. Replace sterile end caps as necessary.
- Reduction in bacterial contamination for antimicrobial connectors has not been shown to correlate with a reduction in infections in patients. Clinical studies to evaluate this have not been performed. Silicone plug contains CAS #7440-22-2 (silver-magnesium-sodium-boron-phosphateglass); polycarbonate cannula contains CAS #265467-11-8 (silver-sodium-hydrogen-zirconium-phosphate).
- Reuse negatively impacts performance / sterility potentially resulting in product failure / contamination.
- Slight discoloration of tubing is a normal result of sterilization and does not affect the safety or sterility of the set.
- Federal (USA) law restricts this device to sale by or on the order of a physician.
- This device should be changed in accordance with current, recognized guidelines of IV therapy.
- When using a set with a MicroClave in the middle of the assembly for high pressure procedures assure that the high pressure infusion is conducted through the MicroClave in the middle of the assembly and the upstream clamp is engaged to isolate other upstream components from high pressures.

COMPONENT	TO PRIME	TO USE	CAUTIONS
Spiros Closed Male Luer	<ul style="list-style-type: none"> • Attach priming cap to Spiros and verify flow. • Once priming of the set and components is complete, remove cap. 	<ul style="list-style-type: none"> • Spiros is normally closed until it is fully engaged with a female luer or needlefree connector. • Attach to mating device by pushing and twisting until secure. 	<ul style="list-style-type: none"> • Do not use luer style end cap on Spiros as this will activate the connector and allow flow.
Clave MicroClave NanoClave Antimicrobial Clave Antimicrobial MicroClave	<ul style="list-style-type: none"> • Prime the connector and set by attaching priming device to connector and inverting to expel air. • Repeat for multiple extensions as necessary. 	<ul style="list-style-type: none"> • Disinfect connector and allow to dry. • Attach syringe or infusion set by pushing straight into the connector and twisting until a friction fit is achieved. Do not overtighten as this could damage the silicone.. • For sets/devices containing a rotating collar, ensure the luer connector is fully engaged prior to tightening the collar. • Flush Clave in accordance with facility protocol. • Flush MicroClave/NanoClave with normal saline or in accordance with facility protocol. • Silicone septum may become displaced or damaged if significant backpressures are applied. • It is not recommended that connectors be changed after use with blood products; devices should be flushed after use with blood in accordance with facility protocol. 	<ul style="list-style-type: none"> • Do not use needles, blunt cannulas or luer caps on connectors. • Connectors are compatible with ISO male luers having an internal diameter between 0.062" and 0.110". Access connectors straight on. (without angled or sliding entry). • MicroClave/NanoClave contains polycarbonate
1o2® Valves	<ul style="list-style-type: none"> • Loosen cap on side port, depress button briefly and observe flow out of port. • Tighten side cap, or remove and attach primed administration device. 	<ul style="list-style-type: none"> • To Infuse: Attach administration device to side port and infuse. Sideport requires pressure to activate. • To Aspirate: attach syringe, depress and hold button while drawing back on syringe plunger. When aspiration is complete, release button and valve automatically returns to closed position. • Recap port with sterile cap after use. 	<ul style="list-style-type: none"> • Do not leave side ports exposed, replace sterile end caps as necessary.
Stopcocks Clave Stopcock NanoClave Stopcock	<ul style="list-style-type: none"> • Use directional turn handle to set-up desired action. • Loosen cap to allow flow. • For stopcocks w/ valves or needlefree connectors on sideports, infuse to prime. 	<ul style="list-style-type: none"> • Use directional turn handle to set-up desired action. • Remove end cap to infuse or aspirate. • When complete turn handle OFF to sideport and replace end cap. 	<ul style="list-style-type: none"> • Do not leave side ports exposed, replace sterile end cap as necessary.
In-line Filters	<ul style="list-style-type: none"> • Orient filter such that inlet is vertically above the outlet. • Filter will automatically prime in the vertical position. 	<ul style="list-style-type: none"> • PN solutions (2-in-1) without lipids can be administered using a 0.22 micron filter. Sets with these filters should be changed at least every 72-hours. • IVFE or TNA (3-in-1) solutions containing lipids must be administered using a 1.2 micron filter. Sets with these filters should be changed at least every 24-hours. • Performance may be impacted by, use of alcohol; nature of fluid used, concentration and duration of use. 	<ul style="list-style-type: none"> • A filter that clogs during infusion should be replaced. Attempting to clear the filter with pressure may lead to breakage.
Eeasydrop® Flow Regulator Eeasydrop double scale flow regulator has one scale for liquids having density lower than 10% (light solutions) and a second scale for liquids having higher density (40% scale).	<ul style="list-style-type: none"> • Easydrop is provided in the OPEN position, verify fluid is passing through. • Turn Easydrop to OFF position and verify fluid has stopped. • If fluid is dripping when in OFF position, replace unit. 	<ul style="list-style-type: none"> • Use of a 15 micron filter on administration set is suggested in order to prevent crystals from blocking the Easydrop. • Adjust Easydrop to the desired rate. • Dial numbers are approximate. Drop counting is necessary to confirm proper flow rate. High viscosity solutions can cause lower flow rates than indicated on the scale; increase head height to compensate. • To change the flow rate adjust height of the solution container. Raise the container to increase flow rate, lower the container to decrease flow rate. • To stop flow, close clamp and turn Easydrop to OFF position. 	<ul style="list-style-type: none"> • Do not use Easydrop in OPEN position, this will cause an uncontrolled delivery (\approx3 liters/hr). • Blood, emulsions, or medications not totally soluble in the carrier solution can not be administered through the Easydrop regulator. • Drop counting is necessary to confirm flow rate. If necessary, adjust the height of the IV solution container to increase or decrease flow rate. • To assure OFF, turn Easydrop to OFF position and close clamp. • Administer bolus injection below controller only.

IN-BX199020460.0

Do Not Reuse

Do Not Resterilize

Non pyrogenic fluid path

Do not use if package is opened or damaged

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Directions for Use

Single patient use, do not resterilize. Sterile, non-pyrogenic fluid path in unopened undamaged package.

DIRECTIONS:

1. Using aseptic technique, open package and tighten connections where applicable.
2. Once primed, attach to patient's catheter.
3. Close clamp and chamber vent cap. If using an Add-On set, insert spike of administration set into spike adapter.
4. Prepare IV solution container in accordance with manufacturer's instructions. Remove end cap on set spike and insert into solution container. (Repeat for multiple spikes if necessary.) For sets containing drip chambers with multiple drip rates, select spike with appropriate drip rate for type of infusion.
5. Hang solution container in accordance with facility protocol.
 - **Drip Chamber:** Squeeze drip chamber to approximately 1/2 full. Fill chamber with ball valve to 3/4 full (~4.73mL). Do not overfill.
 - **Spike Only:** Open clamp and prime.
6. To prime and use various optional components, open clamps on set and refer to chart below.
7. Once primed attach, to desired patient's catheter or compatible y-site on the primary set. For sets with rotating luers, ensure the luer is fully engaged prior to locking down the collar.
8. Use roller clamp(s) to adjust flow rate as desired.

NOTES:

- Pressure infusion not to exceed 400psig.
- Disinfect all needlefree ports using an aggressive circular motion for three (3) seconds.
- Reduction in bacterial contamination for antimicrobial connectors has not been shown to correlate with a reduction in infections in patients. Clinical studies to evaluate this have not been performed. Silicone plug contains CAS #7440-22-2 (silver-magnesium-sodium-boron-phosphateglass); polycarbonate cannula contains CAS #265467-11-8 (silver-sodium-hydrogen-zirconium-phosphate).
- Do not leave open ports/hubs exposed. Replace sterile end caps as necessary.
- Reuse negatively impacts performance / sterility potentially resulting in product failure / contamination.
- Slight discoloration of tubing and/or drip chamber is a normal result of sterilization and does not affect the safety or sterility of the set.
- Federal (USA) law restricts this device to sale by or on the order of a physician.
- This device should be change in accordance with current, recognized guidelines of IV therapy.

	10 Drop ML/Hour	15 Drop ML/Hour	20 Drop ML/Hour	60 Drop ML/Hour
10 mL	1.5 DPM	2.5 DPM	3 DPM	10 DPM
30 mL	5 DPM	7.5 DPM	10 DPM	30 DPM
50 mL	8.5 DPM	12.5 DPM	17 DPM	50 DPM
80 mL	13.5 DPM	20 DPM	27 DPM	80 DPM
100 mL	16.5 DPM	25 DPM	33 DPM	100 DPM
125 mL	21 DPM	31 DPM	42 DPM	125 DPM
150 mL	25 DPM	37.5 DPM	50 DPM	150 DPM
200 mL	33.5 DPM	50 DPM	67 DPM	200 DPM
* Fractions are rounded to closest number except where shown.				

COMPONENT	TO PRIME	TO USE	CAUTIONS
Burette Chamber	<ul style="list-style-type: none"> Close clamp below burette. Ensure air inlet is open. Fill burette with approximately 30 mL fluid. 	<ul style="list-style-type: none"> For continuous infusion, close clamp on set, open burette slide clamp and fill burette chamber to desired level. Close burette air inlet port. Ensure burette slide clamp is open. Adjust flow rate with roller clamp on set. To add fluid or medications to burette, close roller clamp, ensure burette airway inlet port is open. Add medication through burette valve port. Mix thoroughly. 	<ul style="list-style-type: none"> Squeezing burette may cause breakage.
Hand Pump	<ul style="list-style-type: none"> Invert and fill pump chamber completely. 	<ul style="list-style-type: none"> For pressure/bolus infusions open clamp fully. Alternately squeeze and release pump chamber allowing it to fill completely between push action. Discontinue pumping when container is empty and close clamp. 	
Spiros® Closed Male Luer	<ul style="list-style-type: none"> Attach priming cap to Spiros and verify flow. Once priming of the set and components is complete, remove cap. 	<ul style="list-style-type: none"> Spiros is normally closed until it is fully engaged with a female luer or needlefree connector. Attach to mating device by pushing and twisting until secure. 	<ul style="list-style-type: none"> Do not use luer style end cap on Spiros as this will activate the connector and allow flow.
Clave® MicroClave® NanoClave® Antimicrobial Clave® Antimicrobial MicroClave®	<ul style="list-style-type: none"> Prime the connector and set by attaching priming device to connector and inverting to expel air. Repeat for multiple extensions as necessary. 	<ul style="list-style-type: none"> Disinfect connector and attach administration device by pushing straight into the connector and twisting until secure. For sets/devices containing a rotating collar, ensure the luer connector is fully engaged prior to tightening the collar. Flush Clave in accordance with facility protocol. Flush MicroClave/NanoClave with normal saline or in accordance with facility protocol. It is not recommended that connectors be changed after use with blood products; devices should be flushed after use with blood in accordance with facility protocol. 	<ul style="list-style-type: none"> Do not use needles, blunt cannulas or luer caps on connectors. Connectors are compatible with ISO male luers having an internal diameter between 0.062" and 0.110". Access connectors straight on. (without angled or sliding entry). MicroClave/NanoClave contains polycarbonate.
1o2® Valves	<ul style="list-style-type: none"> Loosen cap on side port, depress button briefly and observe flow out of port. Tighten side cap, or remove and attach primed administration device. 	<ul style="list-style-type: none"> To Infuse: Attach administration device to side port and infuse. Sideport requires pressure to activate. To Aspirate: attach syringe, depress and hold button while drawing back on syringe plunger. When aspiration is complete, release button and valve automatically returns to closed position. Recap port with sterile cap after use. 	<ul style="list-style-type: none"> Do not leave side ports exposed, replace sterile end caps as necessary.
Stopcocks Clave Stopcock NanoClave Stopcock	<ul style="list-style-type: none"> Use directional turn handle to set-up desired action. Loosen cap to allow flow. For stopcocks w/ valves on sideports, infuse to prime. 	<ul style="list-style-type: none"> Use directional turn handle to set-up desired action. Remove end cap to infuse or aspirate. When complete turn handle OFF to sideport and replace end cap. 	<ul style="list-style-type: none"> Do not leave side ports exposed, replace sterile end cap as necessary.
In-line Filters Blood Filters	<ul style="list-style-type: none"> Orient filter such that inlet is vertically above the outlet. Filter will automatically prime in the vertical position. Squeeze and release blood filter until filter is covered. 	<ul style="list-style-type: none"> PN solutions (2-in-1) without lipids can be administered using a 0.22 micron filter. Sets with these filters should be changed at least every 72-hours. IVFE or TNA (3-in-1) solutions containing lipids must be administered using a 1.2 micron filter. Sets with these filters should be changed at least every 24-hours. 	<ul style="list-style-type: none"> A filter that clogs during infusion should be replaced. Attempting to clear the filter with pressure may lead to breakage.
Easydrop® Flow Regulator Easydrop double scale flow regulator has one scale for liquids having density lower than 10% (light solutions) and a second scale for liquids having higher density (40% scale).	<ul style="list-style-type: none"> Easydrop is provided in the OPEN position, verify fluid is passing through. Turn Easydrop to OFF position and verify fluid has stopped. If fluid is dripping when in OFF position, replace unit. 	<ul style="list-style-type: none"> Use of a 15 micron filter on administration set is suggested in order to prevent crystals from blocking the Easydrop. Adjust Easydrop to the desired rate. Dial numbers are approximate. Drop counting is necessary to confirm proper flow rate. High viscosity solutions can cause lower flow rates than indicated on the scale; increase head height to compensate. To change flow rate, adjust height of the solution container. Raise the container to increase flow rate, lower the container to decrease flow rate. To stop flow, close clamp and turn Easydrop to OFF position. 	<ul style="list-style-type: none"> Do not use Easydrop in OPEN position, this will cause an uncontrolled delivery (~3 liters/hr). Blood, emulsions, or medications not totally soluble in the carrier solution can not be administered through the Easydrop regulator. Drop counting is necessary to confirm flow rate. If necessary, adjust the height of the IV solution container to increase or decrease flow rate.
Stat2® Flow Controller STAT2 is the registered trademark of Conmed Corporation	<ul style="list-style-type: none"> Turn STAT2 to OPEN position, invert and tap to expel air. Turn STAT2 to OFF position and verify fluid has stopped. If fluid is dripping when in OFF position, replace unit. 	<ul style="list-style-type: none"> Adjust the dial to desired rate. Dial numbers are approximate. Drop counting is necessary to confirm proper flow rate. To stop flow, close clamp and turn STAT2 dial to OFF position. 	<ul style="list-style-type: none"> Blood, emulsions, or medications not totally soluble in the carrier solution can not be administered through the STAT2 Controller. Drop counting is necessary to confirm flow rate. If necessary, adjust the height of the IV solution container to increase or decrease flow rate. To assure OFF, turn STAT2 dial to OFF position and close clamp. Administer bolus injection below controller only.