

Baxter**ONE-LINK****Non-DEHP Standard Bore Catheter Extension Set**

Needle-Free IV Connector with Neutral Fluid Displacement, Power Injectable (325 psi / 2241 kPa)

7.6" (19.3 cm), 0.59 mL

Male Luer Lock Adapter with Retractable Collar

Sterile, Nonpyrogenic.

**Needle-free IV Connector****Do not use if package has been opened or damaged.****Indications for Use:**

The Baxter ONE-LINK Needle-free IV connector is intended for single patient use with a vascular access device for the administration of drugs and solutions without needles, thus eliminating the potential for needle-stick injuries during use. This device is an in-line injection site which can be connected to standard male Luer adapters (e.g., syringe or sets) for continuous or intermittent fluid administration or withdrawal of fluids. This device may be used with low pressure power injectors.

Directions: Use aseptic technique.

1. Prior to power injection, ensure all parts of the I.V. system are compatible. **Tubing and slide clamp are labeled with 325 psi / 2241 kPa text as visual indicators of compatibility with power injection.**
2. **Secure in-line Luer connection.** Do not overtighten. Remove protector as required.
3. When using ONE-LINK connector, grasp finger grip to stabilize.
4. **Swab** Luer activated surface with preferred antiseptic **prior to first use and before every subsequent connection.**
5. Ensure slide clamp is open.
6. Prime set, purge air.
7. For each access, firmly push male Luer of syringe or administration set directly against ONE-LINK connector's Luer activated surface and rotate until connection is secure. Engage Luer locking features, if applicable.
8. Hold the ONE-LINK connector securely while detaching syringe or administration set.
9. Attach male Luer adapter to vascular access device using a firm push and twist motion and then engage the Luer lock collar to prevent accidental disconnection.
10. Open and close slide clamp as required. Clamp extension set when not in use.

Cautions:

- Ensure all in-line Luer connections are secure prior to use.
- To avoid the risks associated with inadvertent ONE-LINK connector disconnection, clamp extension set when not in use.
- Prior to power injection, ensure all parts of the I.V. system are compatible. Replace if a pressure over 325 psi / 2241 kPa is applied to the set.
- Do not allow air to be trapped in set.
- Puncturing set components may cause air embolism.
- Do not regulate flow with slide clamp.
- **Use of Luer lock connection is recommended. If Luer slip connection is used, insert into ONE-LINK connector using a firm push and twist motion. Do not leave Luer slip connection unattended.**
- **Do not use ONE-LINK connector with needles or cannulas. Attempting such access will render the product damaged; replace immediately.**
- If ONE-LINK connector cannot be cleared of blood after blood infusion/sampling, replace immediately. If blood remains on ONE-LINK connector's Luer activated surface after sampling, swab surface.
- Trace lines before connection. Do not connect any compressed gas devices to vascular injection sites.
- Replace the ONE-LINK connector if exposed to gross contamination.
- Replace the ONE-LINK connector if damaged.
- USA: Rx only.
- Single use only. Do not resterilize. Reuse or reprocessing of a single use device may lead to contamination and compromised device function or structural integrity.

Notes:

- a) This product is not made with natural rubber latex.
- b) This product does not contain DEHP.
- c) This product does not contain metal components.
- d) This set can be used to administer lipid-based solutions.
- e) Flush ONE-LINK connector after injection to prevent inadvertent mixing of incompatible medications/fluids.
- f) Flush ONE-LINK connector with 10 mL or more after blood infusion/sampling. If 10 mL flush cannot be performed after blood infusion/sampling, replace ONE-LINK connector immediately.
- g) Luer lock blood collection devices, syringes and accessories are recommended for blood sampling.
- h) Always protect the ONE-LINK connector with an impermeable cover during showering or bathing.
- i) If intermittently disconnecting the administration set from the ONE-LINK connector, immediately cover the male Luer of the administration set with a sterile replacement protector.
- j) Replace per institutional protocol. The ONE-LINK Needle-free IV connector can be used for up to 200 actuations and over a period of 7 days; replace device whichever comes first. The extension set is qualified for use over a period of 7 days; replace after 7 days of use.
- k) Volume and length are approximate.
- l) This device may be used with low pressure power injectors having a maximum pressure of 325 psi / 2241 kPa and a maximum flow rate of 10 mL/s.

For Product Information 1-800-933-0303

Manufactured by an affiliate of
Baxter Healthcare Corporation
 Deerfield, IL 60015 USA
 Made in Costa Rica

Baxter and One-Link are trademarks of Baxter International Inc.

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Baxter et One-Link sont des marques de commerce de Baxter International Inc.