## EpiFuse® Catheter Connector

en English  
EpiFuse® Catheter Connector

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

**NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.**

**1. DESCRIPTION:**  
The Portex® EpiFuse® Catheter Connector is a sterile, single use, reversible snap style connector. It is designed to provide a secure interface between the catheter and the delivery system and is designed to minimize catheter occlusion. The connector is intended for catheters using 16G, 17G and 18G needles.

**2. INDICATIONS FOR USE:**  
The Portex® EpiFuse® Catheter Connector is indicated for use with catheters commonly found in local or regional anesthesia systems.

**3. CONTRAINDICATIONS:**  
There are no known contraindications for use of this device.

**4. PRECAUTIONS:**  
Follow standard infection control procedures as specified by the Centers for Disease Control and Prevention (USA) or local equivalent.

**5. WARNINGS:**  
Fully insert catheter into connector. Failure to do so may result in connector occlusion or catheter disconnection, resulting in inadequate delivery.

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
**6. INSTRUCTIONS FOR USE:**  
After placing catheter and removing needle, insert exposed end of catheter into connector sleeve as far as possible. See figure 1.

**6.1** Snap catheter connector together until closed. An audible click may be heard indicating complete closure. See figure 2.

**6.3** Continue with procedure.

**6.4** To open a closed connector, insert a male Luer slip fitting into the release aperture until the halves of the connector separate. See figure 3.

**6.5** Dispose of contaminated product in a safe manner according to Centers for Disease Control and Prevention (USA) and Federal/State/local regulations (EPA, OSHA) and health care facility guidelines or local equivalent.

EpiFuse, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.

© 2017, 2019 Smiths Medical. All rights reserved.

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**Manufacturer:**  
Smiths Medical ASD, Inc.

6000 Nathan Lane North  
Minneapolis, MN 55442 USA  
Tel: 1 800 258 5361 (US/CA)  
Tel: +1 614 210 7300

**EC REP** European Representative:  
Smiths Medical Czech Republic a.s.  
Olomoucká 306, Hranice 1 - Město,  
753 01 Hranice, Czech Republic  
Tel: +44 (0) 1233 722100  
www.smiths-medical.com

**en** English

### READ COMPLETELY BEFORE USE

THESE INSTRUCTIONS ARE INTENDED FOR USE WITH THE FOLLOWING SMITHS MEDICAL PRODUCT:

## Regional Anesthesia Catheter Securement Device

REF	Device fits Catheters using needle size(s) below:	
Catheter O.D. (mm)	Needle Gauge	
100/399/216	1.1 mm / 1.0 mm	16/17
100/399/218	0.8 mm	18

**1. DESCRIPTION:**  
The security device has an adhesive backing to facilitate skin fixation and a clear central inspection window to allow observation of the insertion site.

**2. INDICATIONS FOR USE:**  
A sterile single use securement device intended for securing regional anesthesia catheters at the skin insertion site to help prevent catheter migration.

### 3. WARNINGS:

3.1 Failure to direct the catheter through the groove on the small central bridge may lead to occlusion of the catheter and failure to deliver medication.

### 4. PRECAUTIONS:

4.1 Ensure that the catheter insertion site is clean and dry before applying the clamp.

4.2 Do not use the securement device in conjunction with aerosol dressings as this may adversely affect adhesion.


4.3 The securement device should be applied using an aseptic technique.

4.4 Avoid applying the securement device in a position where excessive flexing may occur as this may adversely affect adhesion.

4.5 When used to secure epidural catheters, to minimize the risk of catheter migration, particularly in obese patients, at least 4 cm of catheter should be inserted into the epidural space and the patient's back should be in the deflexed position before applying the securement device.

4.6 The maximum period of use is 7 days. (Consult catheter IFU for suggested duration of catheter use. When catheter is removed or changed, follow LOCKIT PLUS™ removal instructions below).

4.7 Blood or bodily fluid may obscure the site of insertion. This can be addressed by opening the "window" and cleaning using standard institutional protocols or by removing the device and replacing it with a new one after cleaning the site.

**5. INSTRUCTIONS FOR USE:**  
5.1 Once the catheter has been successfully positioned before attaching the connector, open the securement device fully and feed the catheter up through the central hole as shown in **Diagram 1**.

5.2 Remove the backing paper and, supporting the catheter, position the securement device onto the skin centrally over the insertion site.

5.3 Press the pad firmly down to ensure good skin contact.

5.4 Locate the catheter into the groove on the small central bridge and push the securement device closed fully until it "clicks" in the latched position (**Diagram 2**). Cover the securement device with a sterile clear waterproof dressing.

- 5.5 **LOCKIT PLUS™ removal:**
    - 5.5.1 Unlatch the device.
    - 5.5.2 Gently remove adhesive layer from patient according to facility protocols for dressing removal.
  - 5.6 After single use, dispose of contaminated product in a safe manner according to Centers for Disease Control and Prevention, USA and/or State/Local regulations (EPA, OSHA) and health care facility guidelines or local equivalent.
- LOCKIT PLUS, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.
- © 2017, 2019 Smiths Medical. All rights reserved.

<b>nl</b>	Gebruiksaanwijzing
<b>sv</b>	Bruksanvisning
<b>da</b>	Brugsanvisning
<b>no</b>	Bruksanvisning
<b>fi</b>	Käytöohjeet
<b>el</b>	Oδηγίες χρήσης

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4.17 Do not use the securement device in conjunction with aerosol dressings as this may adversely affect adhesion.

4.18 The securement device should be applied using an aseptic technique.

4.19 Avoid applying the securement device in a position where excessive flexing may occur as this may adversely affect adhesion.

4.20 When used to secure epidural catheters, to minimize the risk of catheter migration, particularly in obese patients, at least 4 cm of catheter should be inserted into the epidural space and the patient's back should be in the deflexed position before applying the securement device.

4.21 The maximum period of use is 7 days. (Consult catheter IFU for suggested duration of catheter use. When catheter is removed or changed, follow LOCKIT PLUS™ removal instructions below).

4.22 Blood or bodily fluid may obscure the site of insertion. This can be addressed by opening the "window" and cleaning using standard institutional protocols or by removing the device and replacing it with a new one after cleaning the site.

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5. INSTRUCTIONS FOR USE:

5.1 Once the catheter has been successfully positioned before attaching the connector, open the securement device fully and feed the catheter up through the central hole as shown in **Diagram 1**.

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**(latched)**



**5. INSTRUCTIONS FOR USE:**

5.1 Once the catheter has been successfully positioned before attaching the connector, open the securement device fully and feed the catheter up through the central hole as shown in **Diagram 1**.

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5.4 Locate the catheter into the groove on the small central bridge and push the securement device closed fully until it "clicks" in the latched position

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Epifuse, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries.
EpiFuse, znak towarowy Portex oraz znak towarowy Smiths Medical są znakami handlowymi Smiths Medical. Symbol ® oznacza, że znak towarowy jest zarejestrowany w Urzędzie Patentowym Stanów Zjednoczonych (U.S. Patent and Trademark Office) i urzędach patentowych innych krajów.
6.6 Epifuse, znak towarowy Portex oraz znak towarowy Smiths Medical są znakami handlowymi Smiths Medical. Symbol ® oznacza, że znak towarowy jest zarejestrowany w Urzędzie Patentowym Stanów Zjednoczonych (U.S. Patent and Trademark Office) i urzędach patentowych innych krajów.
LOCKIT PLUS, designmerket Portex og designmerket Smiths Medical er varemærker som tilhører Smiths Medical. Symbolet ® angir at varemærket er registrert ved U.S. Patent and Trademark Office og i visse andre land.