CytoRich is an autologous blood-derived product comprising a regenerative component and an anti-inflammatory/catabolic component. It contains biologically active proteins, including PDGF, TGF-beta, IL-1ra, and TIMP-1,2. These Components play a critical role in tissue repair, inflammation modulation, and extracellular matrix regulation. The efficacy of the **CytoRich** procedure depends on proper same-day protocol preparation.

Indications for Use: For autologous use only. **DO NOT** use any component of the system if the packaging is opened or damaged. **DO NOT** clean and/or re-sterilize. Single-use only. CytoRich XL kit is designed for the treatment of tendon and joint inflammation as well as for hair loss therapy.

STEP	PROCEDURE DESCRIPTION
1 Blood draw	Add 3 mL of 4% Sodium Citrate to 60 mL syringe. Before puncture, clean the site with 2% chlorhexidine gluconate in 70% isopropyl alcohol solution, allow drying for 30 seconds. Then slowly draw 57mL of blood, total amount should reach 60mL. Gently mix anticoagulant with blood in syringe, tipping up and down several times. Clean the tops of the tubes with 2% chlorhexidine gluconate in 70% isopropyl alcohol solution; allow to dry for 30 seconds. Attach 18 gauge needle to the syringe and transfer blood to the 6 tubes provided.
2 Preparation of regenerative portion	Place three tubes in the centrifuge and spin using program N1. When the spin is complete, carefully remove the tube caps without touching the edges to prevent contamination of the vial. Using a 10 mL syringe with an 18-gauge needle attached, carefully withdraw the yellow plasma portion, avoiding the collection of sedimented erythrocytes and leukocytes. The syringes should be labeled with patient identifiers and stored at room temperature for up to 6 hours.
3 Preparation of anti-inflammatory/catabolic portion	Store the remaining three untouched tubes at room temperature for 3.5 to 6 hours. Do not exceed 6 hours of holding time. Before product preparation, inspect the color and physical appearance of the contents prior to centrifugation. If there is any indication or suspicion of contamination—such as a green color or cloudy appearance—the product must be discarded immediately. Once the tubes have reached 3.5 to 6 hours of holding time, retrieve them and centrifuge using program N2. When centrifugation is complete, carefully remove the tube caps without touching the edges to prevent contamination of the vial. Using a 10 mL syringe with an 18-gauge needle attached, carefully withdraw the yellow portion avoiding the collection of sedimented erythrocytes and leukocytes. These syringes must be labeled with patient identifiers.
4 Preparation of final CytoRich composition	Attach a sterile filter to the regenerative portion syringe. Connect the sterile connector to the opposite side of the filter. Attach a sterile 10 mL syringe to the connector. Pass up to 5 mL of plasma through the filter and connector into the receiving syringe. This step mechanically activates the regenerative portion. Remove and discard the filter. Set aside the syringe with activated plasma. Retain the syringe with the remaining inactivated plasma, and repeat the above steps using a new sterile filter, connector, and syringe to activate the remaining plasma. To create the final CytoRich product, attach a new sterile filter to the anti-inflammatory/catabolic portion syringe. Pass up to 5 mL of plasma through the filter into the syringe containing the previously activated plasma. This creates the final CytoRich composition with a 1:1 ratio, totaling 10 mL per syringe. Repeat this step for the second syringe containing activated plasma. Discard the used filter and connector. Cap each syringe with a sterile needle. Discard any leftover plasma not used in the final product.

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


Surshield® Safety Winged Infusion and Blood Collection Sets

Terumo® Surshield® Safety Winged Needle Sets provide:

- One-handed*, easy-to-activate safety device
- Minimal steps for activation
- No change in needle insertion technique required
- Confirmation that needle is securely locked
- Minimize the risk of needlestick injuries

*Note: OSHA bloodborne pathogen standard 1910.1030(d)(2)(vii)(B)

Surshield® Safety Winged Needle Sets have the soft, flexible, interlocking wings that readily conform to body's contour. The needles are bonded to the wings at just under a 3° downward angle, for less risk of vessel damage or occlusion of the lumen. The coil-resistant tubing reduces recoil during use. Surshield® Safety Needle sets are DEHP free and have no components made of natural rubber latex.

Surshield® Safety Winged Blood Collection Sets

Packed: 50/box, 10 boxes/case (500)

Product Code	Description
(12" Tubing)	
MN*SVS21B30	21G x ¾" needle
MN*SVS23B30	23G x ¾" needle
MN*SVS25B30	25G x ¾" needle

(7" Tubing)	
MN*SVS21B18	21G x ¾" needle
MN*SVS23B18	23G x ¾" needle

Surshield® Safety Winged Infusion Sets

Packed: 50/box, 10 boxes/case (500)

Product Code	Description
(3½" Tubing)	
SV*S19BLS	19G x ¾" Needle
SV*S21BLS	21G x ¾" Needle
SV*S23BLS	23G x ¾" Needle
SV*S25BLS	25G x ¾" Needle

(12" Tubing)	
SV*S19BL	19G x ¾" Needle
SV*S21BL	21G x ¾" Needle
SV*S23BL	23G x ¾" Needle
SV*S25BL	25G x ¾" Needle

Multi-Sample Luer Adapter

Sterile and disposable, adapters fit into standard infusion sets transforming them into winged blood collection sets. This allows multiple draws into evacuated tubes. Accepts standard tube holders.

Packed: 100/box, 30 boxes/case (3,000)

Product Code	Description
XX*MN2000TE	VENOJECT® Multi-Sample Luer Adapter (male).

Call your Terumo representative at:

1.800.888.3786

or visit www.TerumoTMP.com

and click on Find Your Local Terumo Rep.

Help shield your staff
against accidental needlesticks
with Surshield®

Surshield® Safety Winged Infusion and Blood Collection Sets

Based on IFU Rev. October 2008.

INSTRUCTIONS FOR USE

Aseptic technique, proper skin preparation and continued protection of the site are essential. Observe Universal Precautions on ALL patients.

PRECAUTION: Keep hands behind the needle at all times during use and disposal.

NOTE: Use only one set for I.V. Administration or Blood Collection. Same set should not be used for both procedures.

1. Open the package.

For I.V. administration: Remove luer adapter from the hub of the blood collection set. Connect an administration set, syringe, etc. to hub and secure.

Proceed to step 5.

NOTE PRIMING VOLUME= 0.41ml

For Blood Collection:

2. Prior to removing the luer adapter protector, tighten the connection of luer adapter to the hub of the blood collection set.

3. Remove protector from luer adapter.

4. Thread into holder. (Fig. 1)

5. Flip the safety shield back away from the needle towards the tubing. Grasp wings.

6. Remove the needle protector. **Caution:** Care should be taken not to touch the needle.

7. Perform venipuncture and confirm proper positioning in the vessel, and proceed with blood collection or I.V. administration per facility protocol.

Note: The draw volume in the initial evacuated blood collection tube will be reduced by approximately 0.5 ml due to the deadspace in the collection set tubing.

AFTER USE

8. Remove tape, if present, from wings.

9. Flip the safety shield forward toward needle. Grasp the wing with your thumb on top of the shield and your index finger under the wing. Completely remove the needle from the puncture site and apply digital pressure to the site using a sterile gauze pad. (Fig. 2)

10. While maintaining the same grip on the device (per step 9 above), activate the safety feature by performing either of the following methods:

Slide your thumb forward over the shield while pinching your finger and thumb together (Fig. 3) (or press the safety shield against a hard surface such as a bedside table) until an audible click is heard and the safety feature is visually confirmed to be activated. (Fig. 4)

Recommended (one-handed) method:

(Reference OSHA 1910.1030 (d)(2)(vii)(B))

OR

Alternate (two-handed) method:

While keeping both hands behind the needle, grasp one of the wing tabs with your free hand to stabilize the device. (Fig. 5) (If necessary, reposition your original grip per step 9 being careful not to press against an exposed area of the needle.) Then, slide your thumb forward over shield while pinching your finger and thumb together (Fig. 3) until an audible click is heard and the safety feature is visually confirmed to be activated. (Fig. 6)

Caution: Needle must be carefully and firmly placed into the lock position of the safety shield.

11. Do not attempt to deactivate the safety device by separating the needle from the safety shield.

12. Dispose of used needles and materials following the policies and procedures of your facility, as well as federal and local regulations for "Sharps Disposal."

These products have no components made of natural rubber latex.

RX ONLY. The advertisement is directed to physicians only, and not to consumers. Refer to product labels and packaging insert for complete warnings, precautions, potential complications, and instructions for use. Products may not have regulatory approval in all countries. Please contact your local sales representative if you have questions about the availability of products in your area.


Fig. 1
Thread into holder


Fig. 2
Remove needle from puncture site, apply gauze


Fig. 3
Pinch activation method - audible click


Fig. 4
Hard surface activation method - audible click


Fig. 5
Grasp one wing to stabilize


Fig. 6
Visually confirm activation

TERUMO MEDICAL CORPORATION 265 Davidson Avenue Somerset NJ 08873 U.S.A.

CONFIDENTIAL/CONTROLLED DOCUMENT - UNCONTROLLED when Printed

TERUMO, TERUMO and SURSHIELD are trademarks owned by Terumo Corporation, Tokyo, Japan, and they are registered with the U.S. Patent & Trademark Office. 10/22 PM-06366

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
GB Instructions for use

INDICATIONS

25 mm Syringe Filters with Supor Membrane are designed for rapid filtration, high throughputs, low protein binding, and broad drug compatibility. For drug preparation and syringe bolus drug administration.

DESCRIPTION

Product Code	Pore Size	Housing
PF2000	0.2 µm	green
PF1200	1.2 µm	clear
PF5000	5 µm	purple

SPECIFICATIONS

Priming Volume: <0.9 mL (approx.)
Retention Volume: <0.1 mL with air purge (drug preparation)

CONTRAINDICATIONS

- Blood and other cellular products will not pass through the filter and may cause blockage.
- This device should not be used to filter solutions that are known to be pyrogenic or contaminated with micro-organisms.
- This filter is not designed to eliminate air. Entrained air upstream of the filter may result in reduced flow. This filter is not intended for intrathecal drug administration.

PRECAUTIONS

- Follow instructions carefully.
- Use aseptic technique.
- Do not use if packaging is opened or damaged.
- When attached to needle, dispose of this device in a puncture resistant container.

CAUTIONS:

- Syringes smaller than 10 mL can generate excessive pressure on filter, which may exceed maximum operating pressure.
- As filter removes particulates, syringe plunger may become harder to use as resistance from filter increases. Change filters when resistance becomes excessive. Failure to change filter may cause housing or membrane rupture as maximum operating pressure is exceeded. This will result in particulate contamination of filtrate.
- Do not aspirate and inject using the same filter.

INSTRUCTIONS FOR USE: Use Aseptic Technique.

ASPIRATION

- Peel backing from blister package.
- Holding package in one hand, secure syringe to female luer of filter using a twisting motion. (See Figure 1)
- Attach a sterile needle or access device to male luer of filter, using a twisting motion.
- Sterilize and penetrate access site in accordance with institutional protocols.
- Gently draw back on syringe plunger to filter desired fluid volume into syringe.
- Remove filter from syringe before injection.

INJECTION

- Fill syringe with solution to be filtered.
- Peel backing from blister package.
- Holding package in one hand, secure syringe to female luer of filter, using a twisting motion. (See Figure 1)
- Attach a sterile needle or access device to male luer of filter, using a twisting motion.
- Sterilize and penetrate access site in accordance with institutional protocols.
- Apply gentle pressure to syringe plunger to inject filtered fluid.

Supor is a registered trademark of Pall Corporation.

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B. Braun Medical Inc.
B. Braun Melsungen AG
Pall International Srl