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Patient Information Leaflet

Sodium Chloride Irrigation Solution 0.9% w/v

1. COMPOSITION

Each ampoule contains:
Sodium Chloride Ph. Eur. 0.9% w/v
Water for Injection q.s.

2. PRESENTATION

- LDPE ampoules
- Pack sizes: 20 ampoules x 15 mL and 25 ampoules x 15 mL

3. MECHANISM OF ACTION

- The Irrigation Solution are intended to rinse the wound whilst also reducing the bacterial load. The antimicrobial effect may be considered as ancillary to the primary mechanical action. Sodium is the major cation of extracellular fluid and functions principally in the control of water distribution, fluid and electrolyte balance and osmotic pressure of body fluids. Chloride, the major extracellular anion, closely follows the physiological disposition of sodium cation in maintenance of acid, base balance, isotonicity and electrodynamic characteristic of the cells.

4. HOW TO OPEN THE DEVICE

- The ampoule shall twist open in one or two turns.

5. INTENDED PURPOSE

- The sterile Irrigation solution is used to wash, flush, or rinse parts of the body (e.g., wounds) or medical supplies/equipment (e.g., dressings, catheters, drainage tubes). This product must not be used for injection or during electrosurgical procedures.

6. INDICATION

- This sterile, non-pyrogenic electrolyte solution is indicated for all general irrigation, washing, rinsing.

7. PRECAUTION

Do not use for irrigation that may result in absorption into the blood.

- NOT FOR INJECTION BY USUAL PARENTERAL ROUTES.
- Do not administer unless solution is clear, seal is intact and container is undamaged.
- Keep out of reach of children
- Discard unused portion.
- We have not studied the effect of Sodium Chloride Irrigation Solution 0.9% w/v on pregnant and breast-feeding women. Hence, we do not recommend to use at this moment.

8. WARNINGS

- FOR IRRIGATION ONLY.
- NOT FOR INJECTION.
- Do not warm above 150°F (66°C).
- After opening container, its contents should be used promptly to minimize the possibility of bacterial growth or pyrogen formation.
- Discard unused portion of irrigating solution since it contains no preservatives.

9. SIDE EFFECTS

- This solution may be absorbed into the blood, especially if it used continuously or if it stays inside your body for a long time. This absorption may cause side effects such as too much fluid in the body. Tell your doctor immediately if any of these unlikely but serious side effects occur: signs of infection (e.g., fever), trouble breathing, swelling, unusual weakness, mental/mood changes (e.g., restlessness, irritability). This is not a complete list of possible side effects.

8. ADVERSE REACTIONS

- Should any adverse reaction occur, discontinue the irritant, evaluate the patient, institute appropriate therapeutic countermeasures and save the remainder of the fluid for examination if deemed necessary.

9. DOSAGE AND ADMINISTRATION

- The dose is dependent upon the capacity or surface area of the structure to be irrigated and the nature of the procedure.

10. OVERDOSAGE

- Symptoms of overdose may include: trouble breathing, swelling of the arms/legs.

11. STORAGE & TRANSPORTATION

- Exposure of products to heat should be minimized.
- Avoid excessive heat.
- Do Not Freeze.
- The recommended storage temperature is 2-32°C.
- Do not warm above 150°F (66°C).
- Protect from light. Do not re-sterilize.
- Discard unused portion.

12. SHELF LIFE: 5 years from the date of manufacture.

13. MANUFACTURED BY

SOLA Pharmaceuticals Canada LP
21 Parr Boulevard
Caledon, Ontario, Canada, L7E 4G3

14. GRAPHICAL SYMBOLS FOR MEDICAL DEVICE

St. No.	Description	Symbols
1.	Lot No.	
2.	Date of Manufacture	
3.	Use by	
4.	Consult Pack Insert	
5.	Do not reuse	
6.	Do not Re-sterile	
7.	CE Marking	
8.	Do not Use if Package is damaged	
9.	Sterilized Using Steam sterilization	
10.	Non Pyrogenic	
11.	Keep Dry	
12.	Keep away from sunlight	
13.	Manufacturer	
14.	Catalogue No.	
15.	Caution	
16.	Temperature limit	
17.	Manufacturer Logo	
18.	Manufacturer Details	

Mfg. Lic. No.: MB09738