

Micro-Kleen3



DNTL-BR-MK3-Q223

Disinfectant Cleaner - Wipes

Micro-Kleen3 cleans and disinfects hard, non-porous surfaces such as acrylic, metals (aluminum, brass, copper, stainless steel, chrome), glass, plastics (polycarbonate, polypropylene, polystyrene and polyvinylchloride (PVC)), Plexiglas, and vinyl. Micro-Kleen3 kills hard-to-kill ESKAPE and ESKAPE MDRO organisms with a 1 minute kill time (see second page.)

Micro-Kleen3 is on the [Health Canada List of disinfectants with evidence for use against COVID-19](#) indicating it can be used against SARS-CoV-2, the virus that causes COVID-19.

Features and Benefits

- **Powerful Scrubbing Wipe.** Our disinfectant wipes have a powerful scrubbing ability to quickly and effectively clean surfaces, while still be safe on surfaces.
- **3 Minutes Kill Time.** Kills tested organisms in 30 seconds to 1 minute. Faster kill times mean faster turnovers.
- **Broad Spectrum.** A true broad spectrum disinfectant that is Bactericidal, Fungicidal, Viricidal, and Tuberculocidal.
- **Surface Safe.** Safe on clinical surfaces such as acrylic, aluminum, brass, copper, stainless steel, chrome, plastics, and more. Maintain clear, bright plastics that are typically clouded by high levels of alcohol.
- **Kills ESKAPE and ESKAPE MDRO.** Kills dangerous ESKAPE and ESKAPE MDRO organisms in 1 minute. ESKAPE bacteria are responsible for approximately 2/3 of all hospital acquired infections.
- **Low Alcohol Formula.** Cleans better than high alcohol formulas (55%). Low alcohol formula (20%) won't affix blood and soil to surfaces.



Order Information ▶ Item Code: MKW06-100 - Size: 100 Wipes (6 x 10 in.) - Case: 6 Units/Case

2060 Winston Park Dr. Unit 400, Oakville, ON L4K 0K5 Canada P 905-825-9300 www.chsltd.com

Micro-Kleen3



Disinfectant Cleaner - Wipes

Directions

Please follow the directions stated on the label and IFU when using this product.

Dispensing Directions

To start feed, remove cover and discard seal. From the center of the wipe roll, pull up a wipe corner, twist it into a point and thread it through the hole located on the container cover. Pull through about one inch. Replace cover. Pull out first wipe and tear off at an angle. Remaining wipes feed automatically, ready for the next use. When not in use, keep center cap closed to prevent evaporation. Keep in an upright storing position.

To Disinfect hard, non-porous surfaces/articles

Take a fresh cleaning wipe and use it to thoroughly wet the hard, non-porous surface that has already been pre-cleaned. Allow surface to remain wet for 3 minutes prior to wiping or allow surface to air dry.

Note: Disinfected surfaces which may come in contact with food must be rinsed with potable water after cleaning and/or disinfection and before reuse.

Kills ESKAPE Organisms

1 minute contact time

- [E] Enterococcus faecium
- [S] Staphylococcus aureus
- [K] Klebsiella pneumoniae
- [A] Acintobacter baumannii
- [P] Pseudomonas aeruginosa
- [E] Enterbacter aerogenes

Kills ESKAPE MDRO Organisms

1 minute contact time

- [E] Enterococcus faecium
- [S] Methicillin-Resistant Staphylococcus aureus (MRSA)
- [K] Carbapenum Resistant Klebsiella pneumoniae NDM
- [A] Multi-Drug Resistant Acintobacter baumannii
- [P] Pseudomonas aeruginosa
- [E] Enterbacter aerogenes