

MicroSure

DUWL Water Quality Testing and Waterlines Treatment

MSH2OS-012Water Quality Testing 12 SyringesMSH2OS-03Water Quality Testing 3 Syringes

- Single use, syringe sampling system water quality indicator.
- Detects bacteria, fungi and biofilm for compliance with Canadian and CDC DUWL guidelines.
- Cost Efficient
- Fast, Easy, Accurate Results in one hour or less
- Accurate Results in one hour or less at room temperature.
- Available in 12 or 3 individual single-use disposable syringe sampling systems / package
- No incubator purchase required.

MSS-01 MicroSure Shock Waterline Treatment

- MicroSure Shock is a concentrated liquid designed for the treatment and cleaning of DUWL's.
- Min 10 Max 20 minute treatment
- Highly concentrated
- 1:10 dilution
- Easy to Use
- Convenient
- Safe

<image><list-item><list-item><list-item><list-item><list-item><list-item><list-item><list-item>





400-2060 Winston Park Dr, Oakville, Ontario L6H 5R7 General inquiries: dental@chsltd.com | To place an order: dentalorders@chsltd.com Phone: 905-825-9300 | Toll Free: 1-800-461-1423 | Fax: 1-800-234-1536

www.chsltd.com

©2023 Canadian Hospital Specialties Limited